



The high pressure homogenizers

EMULSIONS, LIPOSOMES, DISPERSIONS, PARTICLE SIZE REDUCTION,
DEAGGLOMERATION, CELL-DISRUPTION, STABILIZATION AND MORE



Let your projects
take off with
PSI homogenizers

Up to 2000 bar
up to 20 L/h
less than 70dBA
small footprint

**Same as a PSI 20
except for the
price**

The EDU program is open to most educational institutions. Please call your local distributors to see if you qualify. Restrictions apply when reselling EDU instruments for the first five years.

FOOTPRINT

WxLxH 650x900x1126 mm

POWER

3x400V 50Hz / 3x460V 60Hz

2.60 kW / 3.05 kW

5.9 A / 5.9 A

WEIGHT

155 kg



www.psi-instruments.com

Distributor

 Qi technologies

via Monte d'Oro 2/A - 00071 Pomezia (RM)
+ 39 06 9105461 - www.qitech.it