

Its compact size and quiet operation do not compromise its capabilities in terms of maximum pressure and flow rate. The PSI-30 is just great for lab work, with batches as small as 100 ml, but is also an excellent small-scale production workhorse, even in continuous process configurations, with a dead volume of 15ml.





FOOTPRINT	W×L×H 650×900×1126 mm
POWER	3×208V 60Hz, 3×380V 60Hz
	3×220V 60Hz, 3×460V 60Hz
	3×380V 50Hz, 3×400V 50Hz
WEIGHT	162 kg

electropolished shell-tube heat exchanger, an easy to use user interface and built-in data acquisition capabilities, but there are many accessories and options to make the PSI-30 the perfect fit for your needs, like GMP compliant record handling, remote control, NIST traceable temperature and pressure sensors.

Contact us

