


# THE HIGH-PRESSURE HOMOGENIZERS



/ CELL-DISRUPTION / DEAGGLOMERATION / DISPERSIONS  
/ EMULSIONS / LIPOSOMES / PARTICLE SIZE REDUCTION  
/ POLYMER MW REDUCTION / STABILIZATION / ...AND MORE



PSI-INSTRUMENTS.COM

ASK ABOUT OUR  
ACADEMIC PROGRAMS



## PSI-EDU

- / up to 2068 bar
- / up to **50Hz** 20 L/h **60Hz** 23 L/h
- / less than 67 dBA
- / small footprint

**FOOTPRINT** W×L×H  
650×900×1126 mm

**POWER**  
3×380V 50Hz, 3×400V 50Hz  
3×220V 60Hz, 3×480V 60Hz

**WEIGHT**  
155 kg

**FLOWRATE**  
**50Hz** 20 L/h **60Hz** 23 L/h

**MAX PRESSURE**  
2068 bar



**LET YOUR PROJECTS TAKE OFF WITH PSI HOMOGENIZERS.**

**Same as a PSI 20** except for the price the **EDU** program is open to most educational institutions.

Please call your local distributors to see if you qualify.

Restrictions apply when reselling **EDU** instruments for the first five years.

**CONTACT US**  
info@psi-instruments.com  
www.psi-instruments.com

**PSI**®  
PSI Instruments

WORK SMARTER